




# SLFN5 mouse mAb

|                           |                                                                                                                                   |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                | YP-mAb-11184                                                                                                                      |
| Isotype                   | IgG                                                                                                                               |
| Reactivity                | Human; Mouse                                                                                                                      |
| Applications              | WB                                                                                                                                |
| Gene Name                 | SLFN5                                                                                                                             |
| Protein Name              | SLFN5                                                                                                                             |
| Immunogen                 | Synthesized peptide derived from human SLFN5 AA range: 707-757                                                                    |
| Specificity               | This antibody detects endogenous levels of SLFN5 at Human/Mouse                                                                   |
| Formulation               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| Source                    | Monoclonal, Mouse,IgG                                                                                                             |
| Purification              | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.              |
| Dilution                  | WB 1:500-1:2000                                                                                                                   |
| Concentration             | 1 mg/ml                                                                                                                           |
| Purity                    | ≥90%                                                                                                                              |
| Storage Stability         | -20°C/1 year                                                                                                                      |
| Synonyms                  |                                                                                                                                   |
| Observed Band             |                                                                                                                                   |
| Cell Pathway              | nucleus,                                                                                                                          |
| Tissue Specificity        |                                                                                                                                   |
| Function                  | function:May have a role in hematopoietic cell differentiation.,similarity:Belongs to the Schlafen family.,                       |
| Background                |                                                                                                                                   |
| matters needing attention | Avoid repeated freezing and thawing!                                                                                              |
| Usage suggestions         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images



Western Blot analysis of various cells using SLFN5 mouse mAb