




FA21A mouse mAb

Catalog No	YP-mAb-11198
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM21A
Protein Name	FA21A
Immunogen	Synthesized peptide derived from human FA21A AA range: 839-889
Specificity	This antibody detects endogenous levels of FA21A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Purity Storage Stability	≥90% -20°C/1 year
Storage Stability	
Storage Stability Synonyms	
Storage Stability Synonyms Observed Band	-20°C/1 year
Storage Stability Synonyms Observed Band Cell Pathway	-20°C/1 year
Storage Stability Synonyms Observed Band Cell Pathway Tissue Specificity	-20°C/1 year Early endosome membrane . Cell membrane .
Storage Stability Synonyms Observed Band Cell Pathway Tissue Specificity Function	-20°C/1 year Early endosome membrane . Cell membrane .
Storage Stability Synonyms Observed Band Cell Pathway Tissue Specificity Function Background matters needing	-20°C/1 year Early endosome membrane . Cell membrane . similarity:Belongs to the FAM21 family.,

Products Images



Western Blot analysis of various cells using FA21A mouse mAb