




ZBED6 mouse mAb

Catalog No	YP-mAb-11208
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZBED6
Protein Name	ZBED6
Immunogen	Synthesized peptide derived from human ZBED6 AA range :130-190
Specificity	This antibody detects endogenous levels of ZBED6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Nucleus, nucleolus . Cytoplasm . Located predominantly in the nucleolus but is also dispersed throughout the nucleus (By similarity). Mainly cytoplasmic in insulin- and glucagon-positive islet pancreatic cells, however nuclear localization is evidence in some pancreatic endocrine cells (PubMed:24043816). .
Tissue Specificity	Expressed in pancreatic islet cells (at protein level).
Function	
Background	The protein encoded by this transposon-derived intronless gene is a transcriptional repressor that binds to the consensus sequence 5'-GCTCGC-3'. The encoded protein has been shown to repress IGF2 transcription. This gene is located within the first intron of the ZC3H11A gene. [provided by RefSeq, Jul 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZBED6 mouse mAb