




NO66 mouse mAb

Catalog No	YP-mAb-11227
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	NO66 C14orf169 MAPJD
Protein Name	NO66
Immunogen	Synthesized peptide derived from human NO66 AA range: 533-583
Specificity	This antibody detects endogenous levels of NO66 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleolus. Nucleus, nucleoplasm. Granular part of nucleoli. Nucleoplasm, nucleoplasmic foci, some of them associated with nucleoli.
Tissue Specificity	Widely expressed. Overexpressed in lung carcinomas.
Function	function:May play a role in ribosome biogenesis and in the replication or remodeling of certain heterochromatic region.,similarity:Belongs to the MINA53/NO66 family.,similarity:Contains 1 JmjC domain.,subcellular location:Granular part of nucleoli. Nucleoplasm, nucleoplasmic foci, some of them associated with nucleoli.,tissue specificity:Ubiquitous.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NO66 mouse mAb