



C1QT6 mouse mAb

Catalog No	YP-mAb-11338
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C1QTNF6 CTRP6 UNQ581/PRO1151
Protein Name	C1QT6
Immunogen	Synthesized peptide derived from human C1QT6 AA range: 31-81
Specificity	This antibody detects endogenous levels of C1QT6 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using C1QT6 mouse mAb

