




TT39C mouse mAb

Catalog No	YP-mAb-11374
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TTC39C C18orf17
Protein Name	TT39C
Immunogen	Synthesized peptide derived from human TT39C AA range: 229-279
Specificity	This antibody detects endogenous levels of TT39C at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the TTC39 family.,similarity:Contains 3 TPR repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TT39C mouse mAb