




ZBT25 mouse mAb

Catalog No	YP-mAb-11376
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZBTB25 C14orf51 KUP ZNF46
Protein Name	ZBT25
Immunogen	Synthesized peptide derived from human ZBT25 AA range: 249-299
Specificity	This antibody detects endogenous levels of ZBT25 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus.
Tissue Specificity	Expressed mainly in hematopoietic cells and testis.
Function	function:May be involved in transcriptional regulation.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,tissue specificity:Expressed mainly in hematopoietic cells and testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZBT25 mouse mAb