




ZSWM4 mouse mAb

Catalog No	YP-mAb-11423
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZSWIM4
Protein Name	ZSWM4
Immunogen	Synthesized peptide derived from human ZSWM4 AA range: 94-144
Specificity	This antibody detects endogenous levels of ZSWM4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 SWIM-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZSWM4 mouse mAb