




RLF mouse mAb

Catalog No	YP-mAb-11485
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RLF
Protein Name	RLF
Immunogen	Synthesized peptide derived from human RLF AA range: 1397-1447
Specificity	This antibody detects endogenous levels of RLF at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed in fetal and adult tissues.
Function	disease:In some small cell lung carcinoma (SCLC) cell lines, there is an intrachromosomal rearrangements at 1p32 fusing the first exon of the RLF gene with L-myc.,function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 14 C2H2-type zinc fingers.,subunit:Interacts with RIT1 and RIT2.,tissue specificity:Widely expressed in fetal and adult tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RLF mouse mAb