




MISSL mouse mAb

Catalog No	YP-mAb-11559
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	MAPK1IP1L C14orf32
Protein Name	MISSL
Immunogen	Synthesized peptide derived from human MISSL AA range: 150-200
Specificity	This antibody detects endogenous levels of MISSL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the MISS family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MISSL mouse mAb