




ISOC1 mouse mAb

Catalog No	YP-mAb-11585
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ISOC1 CGI-111
Protein Name	ISOC1
Immunogen	Synthesized peptide derived from human ISOC1 AA range: 193-243
Specificity	This antibody detects endogenous levels of ISOC1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	cytoplasm,peroxisome,extracellular exosome,
Tissue Specificity	
Function	similarity:Belongs to the isochorismatase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ISOC1 mouse mAb