




CB016 mouse mAb

Catalog No	YP-mAb-11606
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C2orf16
Protein Name	CB016
Immunogen	Synthesized peptide derived from human CB016 AA range: 1425-1475
Specificity	This antibody detects endogenous levels of CB016 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	nucleus,extracellular exosome,
Tissue Specificity	
Function	domain:The P-S-E-R-S-H-H-S repeats give rise to an antiparallel beta-structure.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CB016 mouse mAb