



EFHD2 mouse mAb

Catalog No	YP-mAb-11611
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	EFHD2 SWS1
Protein Name	EFHD2
Immunogen	Synthesized peptide derived from human EFHD2 AA range: 13-63
Specificity	This antibody detects endogenous levels of EFHD2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane raft . In a mouse immature B-cell line WEHI-231. .
Tissue Specificity	Found in lymphocytes; preferentially expressed in CD8+ cells.
Function	similarity:Contains 2 EF-hand domains.,tissue specificity:Found in lymphocytes; preferentially expressed in CD8+ cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using EFHD2 mouse mAb

