




PCYXL mouse mAb

Catalog No	YP-mAb-11624
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PCYOX1L PSEC0105
Protein Name	PCYXL
Immunogen	Synthesized peptide derived from human PCYXL AA range: 406-456
Specificity	This antibody detects endogenous levels of PCYXL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	
Function	cofactor:FAD.,function:Probable oxidoreductase.,similarity:Belongs to the prenylcysteine oxidase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PCYXL mouse mAb