




PHLP2 mouse mAb

Catalog No	YP-mAb-11626
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PHLPP2 KIAA0931 PHLPP1
Protein Name	PHLP2
Immunogen	Synthesized peptide derived from human PHLP2 AA range: 381-431
Specificity	This antibody detects endogenous levels of PHLP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm. Membrane; Peripheral membrane protein. Nucleus. In colorectal cancer tissue, expression is concentrated in the cytoplasm and nucleus.
Tissue Specificity	In colorectal cancer tissue, expression is highest in the surface epithelium of normal colonic mucosa adjacent to the cancer tissue but is largely excluded from the crypt bases. Expression is lost or significantly decreased in 80% of tested tumors (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PHLP2 mouse mAb