




RBFA mouse mAb

Catalog No	YP-mAb-11732
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RBFA C18orf22
Protein Name	RBFA
Immunogen	Synthesized peptide derived from human RBFA AA range: 286-336
Specificity	This antibody detects endogenous levels of RBFA at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	similarity:Belongs to the rbfA family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RBFA mouse mAb