




K1551 mouse mAb

Catalog No	YP-mAb-11760
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	KIAA1551 C12orf35
Protein Name	K1551
Immunogen	Synthesized peptide derived from human K1551 AA range: 262-312
Specificity	This antibody detects endogenous levels of K1551 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Localizes around gamma-tubulin during M phase. .
Tissue Specificity	
Function	sequence caution:Contaminating sequence.,sequence caution:Contaminating sequence. Potential poly-A sequence.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using K1551 mouse mAb