




P33MX mouse mAb

Catalog No	YP-mAb-11761
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	KIAA1191 P33MONOX
Protein Name	P33MX
Immunogen	Synthesized peptide derived from human P33MX AA range: 186-236
Specificity	This antibody detects endogenous levels of P33MX at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Down-regulated in the occipital lobe of an early stage Alzheimer disease patients.
Function	similarity:Belongs to the UPF0498 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using P33MX mouse mAb