




WSCD1 mouse mAb

Catalog No	YP-mAb-11796
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	WSCD1 KIAA0523
Protein Name	WSCD1
Immunogen	Synthesized peptide derived from human WSCD1 AA range: 442-492
Specificity	This antibody detects endogenous levels of WSCD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the WSCD family.,similarity:Contains 2 WSC domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using WSCD1 mouse mAb