




PAMR1 mouse mAb

Catalog No	YP-mAb-11797
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PAMR1 RAMP FP938 UNQ699/PRO1344
Protein Name	PAMR1
Immunogen	Synthesized peptide derived from human PAMR1 AA range: 57-107
Specificity	This antibody detects endogenous levels of PAMR1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	79kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PAMR1 mouse mAb