




AFAF mouse mAb

Catalog No	YP-mAb-11800
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	AFAF C9orf11
Protein Name	AFAF
Immunogen	Synthesized peptide derived from human AFAF AA range: 223-273
Specificity	This antibody detects endogenous levels of AFAF at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome membrane ; Single-pass type I membrane protein . Cytoplasmic vesicle, secretory vesicle, acrosome inner membrane ; Single-pass type I membrane protein . Cytoplasmic vesicle, secretory vesicle, acrosome outer membrane ; Single-pass type I membrane protein . In the anterior acrosome region, enriched on the inner acrosomal membrane but minimal on the outer acrosomal membrane; in contrast in the posterior acrosome region enriched on both the inner and outer acrosomal membranes. .
Tissue Specificity	Isoform 1 is highly expressed in testis. Isoform 2 is expressed at low levels in skin and blood.
Function	function:Involved in acrosome biogenesis.,tissue specificity:Isoform 1 highly expressed in testis. Isoform 2 is expressed at low levels in skin and blood.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using AFAF mouse mAb