



TEX22 mouse mAb

Catalog No	YP-mAb-11824
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TEX22
Protein Name	TEX22
Immunogen	Synthesized peptide derived from human TEX22 AA range: 5-55
Specificity	This antibody detects endogenous levels of TEX22 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Cytoplasmic vesicle, secretory vesicle, acrosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using TEX22 mouse mAb

