




CA170 mouse mAb

Catalog No	YP-mAb-11826
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C1orf170
Protein Name	CA170
Immunogen	Synthesized peptide derived from human CA170 AA range: 434-484
Specificity	This antibody detects endogenous levels of CA170 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Nucleus . Shows a nuclear localization in the presence of PPARGC1A. .
Tissue Specificity	Muscle-specific expression is increased by endurance exercise.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CA170 mouse mAb