




PTCD3 mouse mAb

Catalog No	YP-mAb-11850
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PTCD3 TRG15
Protein Name	PTCD3
Immunogen	Synthesized peptide derived from human PTCD3 AA range: 624-674
Specificity	This antibody detects endogenous levels of PTCD3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Mitochondrion .
Tissue Specificity	Abundant in testes, skeletal muscle and heart tissue.
Function	similarity:Belongs to the PTCD3 family.,similarity:Contains 10 PPR (pentatricopeptide) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PTCD3 mouse mAb