



LYSM2 mouse mAb

Catalog No	YP-mAb-11857
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LYSMD2
Protein Name	LYSM2
Immunogen	Synthesized peptide derived from human LYSM2 AA range: 100-150
Specificity	This antibody detects endogenous levels of LYSM2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 LysM repeat.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using LYSM2 mouse mAb

