




ZBTB5 mouse mAb

Catalog No	YP-mAb-11874
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZBTB5 KIAA0354
Protein Name	ZBTB5
Immunogen	Synthesized peptide derived from human ZBTB5 AA range: 54-104
Specificity	This antibody detects endogenous levels of ZBTB5 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May be involved in transcriptional regulation.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZBTB5 mouse mAb