




CB083 mouse mAb

Catalog No	YP-mAb-11875
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C2orf83
Protein Name	CB083
Immunogen	Synthesized peptide derived from human CB083 AA range: 56-106
Specificity	This antibody detects endogenous levels of CB083 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	integral component of plasma membrane,
Tissue Specificity	
Function	caution:Although related to reduced folate carrier family, it lacks transmembrane domains, suggesting that it has probably no transporter function.,sequence caution:Translated as Trp.,similarity:Belongs to the reduced folate carrier (RFC) transporter (TC 2.A.48) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CB083 mouse mAb