




TRM13 mouse mAb

Catalog No	YP-mAb-11883
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TRMT13 CCDC76
Protein Name	TRM13
Immunogen	Synthesized peptide derived from human TRM13 AA range: 331-381
Specificity	This antibody detects endogenous levels of TRM13 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	catalytic activity:S-adenosyl-L-methionine + rRNA = S-adenosyl-L-homocysteine + rRNA containing 2'-O-methyluridine.,function:Specifically methylates guanosine-4 in various tRNAs with a Gly(CCG), His or Pro signature.,similarity:Belongs to the TRM13 family.,similarity:Contains 1 CHHC-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TRM13 mouse mAb