




ZCHC7 mouse mAb

Catalog No	YP-mAb-11915
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ZCCHC7 HSPC086
Protein Name	ZCHC7
Immunogen	Synthesized peptide derived from human ZCHC7 AA range: 4-54
Specificity	This antibody detects endogenous levels of ZCHC7 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-2 is the initiator.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 4 CCHC-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZCHC7 mouse mAb