




RTKN2 mouse mAb

Catalog No	YP-mAb-11944
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RTKN2 PLEKHK1
Protein Name	RTKN2
Immunogen	Synthesized peptide derived from human RTKN2 AA range: 473-523
Specificity	This antibody detects endogenous levels of RTKN2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Expressed in lymphocytes, CD4 positive T-cells and bone marrow-derived cells. Also expressed in lung, colon, thymus and brain.
Function	function:May play an important role in lymphopoiesis.,induction:Up-regulated in cells resistant to 25-hydroxy cholesterol (25-OHC). Down-regulated in lymphocytes activated by treatment with phytohemagglutinin (PHA).,similarity:Contains 1 PH domain.,tissue specificity:Expressed in lymphocytes, CD4 positive T-cells and bone marrow-derived cells. Also expressed in lung, colon, thymus and brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RTKN2 mouse mAb