




ISM2 mouse mAb

Catalog No	YP-mAb-11967
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ISM2 TAIL1 THSD3 PSEC0137
Protein Name	ISM2
Immunogen	Synthesized peptide derived from human ISM2 AA range: 65-115
Specificity	This antibody detects endogenous levels of ISM2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	Expressed at high levels in the placenta and at moderate levels in the pancreas, kidney, heart, liver, lung, brain and skeletal muscle.
Function	similarity:Belongs to the isthmin family.,similarity:Contains 1 AMOP domain.,similarity:Contains 1 TSP type-1 domain.,tissue specificity:Expressed at high levels in the placenta and at moderate levels in the pancreas, kidney, heart, liver, lung, brain and skeletal muscle.,
Background	The protein encoded by this gene contains a type 1 thrombospondin domain, which is present in thrombospondin, a number of proteins involved in the complement pathway, as well as in extracellular matrix proteins. Two alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Mar 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ISM2 mouse mAb