




## KDEL2 mouse mAb

|                           |                                                                                                                                   |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                | YP-mAb-11971                                                                                                                      |
| Isotype                   | IgG                                                                                                                               |
| Reactivity                | Human;Rat                                                                                                                         |
| Applications              | WB                                                                                                                                |
| Gene Name                 | KDELC2 UNQ1904/PRO4350                                                                                                            |
| Protein Name              | KDEL2                                                                                                                             |
| Immunogen                 | Synthesized peptide derived from human KDEL2 AA range: 414-464                                                                    |
| Specificity               | This antibody detects endogenous levels of KDEL2 at Human/Rat                                                                     |
| Formulation               | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| Source                    | Monoclonal, Mouse,IgG                                                                                                             |
| Purification              | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.              |
| Dilution                  | WB 1:500-1:2000                                                                                                                   |
| Concentration             | 1 mg/ml                                                                                                                           |
| Purity                    | ≥90%                                                                                                                              |
| Storage Stability         | -20°C/1 year                                                                                                                      |
| Synonyms                  |                                                                                                                                   |
| Observed Band             |                                                                                                                                   |
| Cell Pathway              | Endoplasmic reticulum lumen .                                                                                                     |
| Tissue Specificity        |                                                                                                                                   |
| Function                  | similarity:Belongs to the KDELC family.,similarity:Contains 1 filamin repeat.,                                                    |
| Background                |                                                                                                                                   |
| matters needing attention | Avoid repeated freezing and thawing!                                                                                              |
| Usage suggestions         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images



Western Blot analysis of various cells using KDEL2 mouse mAb