




NOP9 mouse mAb

Catalog No	YP-mAb-11983
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NOP9 C14orf21 KIAA2021
Protein Name	NOP9
Immunogen	Synthesized peptide derived from human NOP9 AA range: 486-536
Specificity	This antibody detects endogenous levels of NOP9 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 6 pumilio repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NOP9 mouse mAb