



# IAH1 mouse mAb

<b>Catalog No</b>	YP-mAb-11984
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human; Mouse;Rat
<b>Applications</b>	WB
<b>Gene Name</b>	IAH1
<b>Protein Name</b>	IAH1
<b>Immunogen</b>	Synthesized peptide derived from human IAH1 AA range: 22-72
<b>Specificity</b>	This antibody detects endogenous levels of IAH1 at Human/Mouse/Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	
<b>Cell Pathway</b>	extracellular exosome,
<b>Tissue Specificity</b>	
<b>Function</b>	function:Probable lipase.,similarity:Belongs to the 'GDGL' lipolytic enzyme family. IAH1 subfamily.,
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

Western Blot analysis of various cells using IAH1 mouse mAb

