




LY10L mouse mAb

Catalog No	YP-mAb-11991
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SP140L
Protein Name	LY10L
Immunogen	Synthesized peptide derived from human LY10L AA range: 178-228
Specificity	This antibody detects endogenous levels of LY10L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	nucleus,
Tissue Specificity	
Function	similarity:Contains 1 HSR domain.,similarity:Contains 1 PHD-type zinc finger.,similarity:Contains 1 SAND domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LY10L mouse mAb