




RERGL mouse mAb

Catalog No	YP-mAb-12014
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RERGL
Protein Name	RERGL
Immunogen	Synthesized peptide derived from human RERGL AA range: 16-66
Specificity	This antibody detects endogenous levels of RERGL at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	intracellular,membrane,
Tissue Specificity	
Function	function:Binds GDP/GTP and may possess intrinsic GTPase activity.,similarity:Belongs to the small GTPase superfamily. Ras family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RERGL mouse mAb