




LRC20 mouse mAb

Catalog No	YP-mAb-12044
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LRRC20 UNQ2429/PRO4989
Protein Name	LRC20
Immunogen	Synthesized peptide derived from human LRC20 AA range: 6-56
Specificity	This antibody detects endogenous levels of LRC20 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	sequence caution:Translated as stop.,similarity:Contains 6 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LRC20 mouse mAb