




ELMD1 mouse mAb

Catalog No	YP-mAb-12108
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ELMOD1
Protein Name	ELMD1
Immunogen	Synthesized peptide derived from human ELMD1 AA range: 73-123
Specificity	This antibody detects endogenous levels of ELMD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	nucleoplasm,cytoplasm,
Tissue Specificity	
Function	similarity:Contains 1 ELMO domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ELMD1 mouse mAb