




WDR53 mouse mAb

Catalog No	YP-mAb-12139
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	WDR53
Protein Name	WDR53
Immunogen	Synthesized peptide derived from human WDR53 AA range: 19-69
Specificity	This antibody detects endogenous levels of WDR53 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the WD repeat WDR53 family.,similarity:Contains 5 WD repeats.,
Background	This gene encodes a protein containing WD domains. The function of this gene is unknown. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using WDR53 mouse mAb