




RNPL1 mouse mAb

Catalog No	YP-mAb-12159
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RNPEPL1
Protein Name	RNPL1
Immunogen	Synthesized peptide derived from human RNPL1 AA range: 112-162
Specificity	This antibody detects endogenous levels of RNPL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	cytoplasm,plasma membrane,
Tissue Specificity	Ubiquitously expressed. Expressed at relatively higher levels in heart and skeletal muscle.
Function	cofactor: Binds 1 zinc ion per subunit.,similarity: Belongs to the peptidase M1 family.,tissue specificity: Ubiquitous. Expressed at relatively higher levels in heart and skeletal muscle.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RNPL1 mouse mAb