



F190B mouse mAb

Catalog No	YP-mAb-12229
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM190B KIAA1128 NPD012
Protein Name	F190B
Immunogen	Synthesized peptide derived from human F190B AA range: 609-659
Specificity	This antibody detects endogenous levels of F190B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton . Associates with microtubules in interphase. Has diffuse expression throughout the cell during mitosis. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using F190B mouse mAb

