



ZFN2A mouse mAb

Catalog No	YP-mAb-12239
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ZFAND2A
Protein Name	ZFN2A
Immunogen	Synthesized peptide derived from human ZFN2A AA range: 42-92
Specificity	This antibody detects endogenous levels of ZFN2A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	similarity:Contains 2 AN1-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using ZFN2A mouse mAb

