




LRRN2 mouse mAb

Catalog No	YP-mAb-12244
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LRRN2 GAC1 LRRN5 UNQ256/PRO293
Protein Name	LRRN2
Immunogen	Synthesized peptide derived from human LRRN2 AA range: 296-346
Specificity	This antibody detects endogenous levels of LRRN2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Overamplified in malignant gliomas.
Function	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.;similarity:Contains 11 LRR (leucine-rich) repeats.;tissue specificity:Overamplified in malignant gliomas.,
Background	The protein encoded by this gene belongs to the leucine-rich repeat superfamily. This gene was found to be amplified and overexpressed in malignant gliomas. The encoded protein has homology with other proteins that function as cell-adhesion molecules or as signal transduction receptors and is a candidate for the target gene in the 1q32.1 amplicon in malignant gliomas. Two alternatively spliced transcript variants encoding the same protein have been described for this gene. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using LRRN2 mouse mAb