




RLA0L mouse mAb

Catalog No	YP-mAb-12262
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RPLP0P6
Protein Name	RLA0L
Immunogen	Synthesized peptide derived from human RLA0L AA range: 20-70
Specificity	This antibody detects endogenous levels of RLA0L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	cytosolic large ribosomal subunit,preribosome, large subunit precursor,
Tissue Specificity	
Function	function:Ribosomal protein P0 is the functional equivalent of E.coli protein L10.,similarity:Belongs to the ribosomal protein L10P family.,subunit:P0 forms a pentameric complex by interaction with dimers of P1 and P2.,subunit:P0 forms a pentameric complex by interaction with dimers of P1 and P2. Interacts with Lassa virus Z protein.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RLA0L mouse mAb