




FA46D mouse mAb

Catalog No	YP-mAb-12267
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM46D
Protein Name	FA46D
Immunogen	Synthesized peptide derived from human FA46D AA range: 84-134
Specificity	This antibody detects endogenous levels of FA46D at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 : 500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	43kD
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	restricted to testis.
Function	similarity:Belongs to the FAM46 family.,
Background	Antibodies against the protein encoded by this gene were found only in plasma from cancer patients. While it may be a target for immunotherapy, the function of this gene is unknown. [provided by RefSeq, Dec 2009],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4°over night