




RBM3 mouse mAb

Catalog No	YP-mAb-12406
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	RBM3 RNPL
Protein Name	RBM3
Immunogen	Synthesized peptide derived from human RBM3 AA range: 69-119
Specificity	This antibody detects endogenous levels of RBM3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cytoplasm . Cell projection, dendrite . Localizes in mRNA granules in dendrites. .
Tissue Specificity	
Function	PTM:Arg-105 is dimethylated, probably to asymmetric dimethylarginine.,similarity:Contains 1 RRM (RNA recognition motif) domain.,
Background	This gene is a member of the glycine-rich RNA-binding protein family and encodes a protein with one RNA recognition motif (RRM) domain. Expression of this gene is induced by cold shock and low oxygen tension. A pseudogene exists on chromosome 1. Multiple alternatively spliced transcript variants that are predicted to encode different isoforms have been characterized although some of these variants fit nonsense-mediated decay (NMD) criteria. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RBM3 mouse mAb