




ALS2CR13 Monoclonal Antibody

Catalog No	YP-mAb-12680
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM117B
Protein Name	Protein FAM117B
Immunogen	The antiserum was produced against synthesized peptide derived from human ALS2CR13. AA range:411-460
Specificity	ALS2CR13 Monoclonal Antibody detects endogenous levels of ALS2CR13 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM117B; ALS2CR13; Protein FAM117B; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein
Observed Band	60kD
Cell Pathway	cytoplasm,
Tissue Specificity	Human salivary gland,Kidney,Skin,
Function	miscellaneous:ALS2CR13 is mapped in the genomic region covering the complete candidate region for Amyotrophic lateral sclerosis 2 (ALS2).,sequence caution:Contaminating sequence. Sequence of unknown origin in the N-terminal part.,
Background	miscellaneous:ALS2CR13 is mapped in the genomic region covering the complete candidate region for Amyotrophic lateral sclerosis 2 (ALS2).,sequence caution:Contaminating sequence. Sequence of unknown origin in the N-terminal part.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ALS2CR13 Monoclonal Antibody