



TMEM101 Monoclonal Antibody

Catalog No	YP-mAb-13703
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMEM101
Protein Name	Transmembrane protein 101
Immunogen	The antiserum was produced against synthesized peptide derived from human TMEM101. AA range:101-150
Specificity	TMEM101 Monoclonal Antibody detects endogenous levels of TMEM101 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMEM101; Transmembrane protein 101; Putative NF-kappa-B-activating protein 130
Observed Band	28kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Lung,Muscle,
Function	function:May activate NF-kappa-B signaling pathways.,
Background	function:May activate NF-kappa-B signaling pathways.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

