




SH2D5 Monoclonal Antibody

Catalog No	YP-mAb-13961
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SH2D5
Protein Name	SH2 domain-containing protein 5
Immunogen	The antiserum was produced against synthesized peptide derived from human SH2D5. AA range:281-330
Specificity	SH2D5 Monoclonal Antibody detects endogenous levels of SH2D5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SH2D5; SH2 domain-containing protein 5
Observed Band	35kD
Cell Pathway	Cell junction, synapse, postsynaptic density .
Tissue Specificity	Hippocampus,Thalamus,
Function	similarity:Contains 1 SH2 domain.,
Background	similarity:Contains 1 SH2 domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SH2D5 Monoclonal Antibody