




Tetraspanin-8 Monoclonal Antibody

Catalog No	YP-mAb-13966
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TSPAN8
Protein Name	Tetraspanin-8
Immunogen	The antiserum was produced against synthesized peptide derived from human TSPAN8. AA range:96-145
Specificity	Tetraspanin-8 Monoclonal Antibody detects endogenous levels of Tetraspanin-8 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TSPAN8; TM4SF3; Tetraspanin-8; Tspan-8; Transmembrane 4 superfamily member 3; Tumor-associated antigen CO-029
Observed Band	26kD
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Gastric, colon, rectal, and pancreatic carcinomas.
Function	similarity:Belongs to the tetraspanin (TM4SF) family.,tissue specificity:Gastric, colon, rectal, and pancreatic carcinomas.,
Background	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins. This gene is expressed in different carcinomas. The use of alternate polyadenylation sites has been found for this gene. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Tetraspanin-8 Monoclonal Antibody